

EASYDRY II

DRYING - WITHIN 7 MINUTES



Device for drying analysis membranes for the determination of gravimetry of Technical Cleanliness of functionally relevant components in accordance with VDA 19.1 and ISO 16232. With significant time and place saving. Minimization of hazards when handling solvents in relation to conventional membrane preparation.

EASYDRY II

SOME HIGHLIGHTS




LED for easy differentiation of six modes with various colors and additional beep e.g. "drying process completed".



REDUCTION of cross-contamination during handling the membranes.




ERGONOMIC WORK TOP POSITION with shelves for an easy and safe membrane handling.



INCREASE OF PROCESS SAFETY during membrane drying by reduction of number of process steps.



EASY MAINTENANCE and reduction of excessive moisture by a slightly emptiable collecting container.



EFFICIENT LABORATORY WORK by double saving of time each at the requested preparation and post-processing of membranes. Drying within 7-9 minutes.

GLÄSER GMBH

Robert-Bosch-Straße 32
D - 72160 Horb am Neckar

Tel.: +49 (0) 7451 53920-0
Fax +49 (0) 7451 53920-44
www.glaeser-group.com
info@glaeser-group.com

